

Product datasheet

Rabbit IgG, polyclonal - Isotype Control ab27472

20 References

Overview

Product name	Rabbit IgG, polyclonal - Isotype Control
Specificity	ab27472 does not react with any antigen of human or animal cells.
Tested applications	Suitable for: IHC-P, IP, Sandwich ELISA
Species reactivity	ab27472 does not react with any species.
General notes	The concentration of this prediluted antibody is 2 µg/ml.

Properties

Form	Prediluted
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.1% Sodium Azide Constituents: PBS, BSA
Purity	Protein G purified
Purification notes	ab27472 is an affinity purified IgG. It was purified from serum collected from a rabbit prior to immunization.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab27472** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent dilution. Use neat for 10 minutes at room temperature. Pretreatment: same as primary antibody.
IP		Use at an assay dependent dilution. PubMed: 21088000
Sandwich ELISA		Use at an assay dependent dilution. PubMed: 19243500

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors