


Product datasheet

Anti-DMRT8.2 antibody ab36912

Overview

<b>Product name</b>	Anti-DMRT8.2 antibody
<b>Description</b>	Rabbit polyclonal to DMRT8.2
<b>Host species</b>	Rabbit
<b>Species reactivity</b>	<b>Predicted to work with:</b> Mouse 
<b>Immunogen</b>	Synthetic peptide conjugated to KLH derived from within residues 1 - 100 of Mouse DMRT8.2. Read Abcam's proprietary immunogen policy (Peptide available as <a href="#">ab36911</a> .)

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Preservative: 0.02% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

<b>Application notes</b>	This antibody gave a positive result in ELISA against the immunizing peptide ( <a href="#">ab36911</a> ). Not yet tested in other applications.
--------------------------	--

Target

<b>Relevance</b>	The vertebrate DMRT gene family encodes putative transcription factors related to the sexual regulators Doublesex (Drosophila) and MAB-3 (Worm). DMRT8 is related to DMRT7 and located on the X chromosome in placental mammals. DMRT8.2 expression is restricted to male and female gonads and may be involved in mouse sexual development.
<b>Cellular localization</b>	Nuclear

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors