

Product datasheet

Anti-Integrin alpha 2b + beta 3 antibody [IVA30] ab3919

6 References

Overview

Product name	Anti-Integrin alpha 2b + beta 3 antibody [IVA30]
Description	Mouse monoclonal [IVA30] to Integrin alpha 2b + beta 3
Host species	Mouse
Tested applications	Suitable for: ELISA, Flow Cyt, IP Unsuitable for: WB
Species reactivity	Reacts with: Cow
Immunogen	Bovine trombocytes.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 15mM Sodium Azide Constituents: PBS, pH 7.4
Purity	>90% by SDS-PAGE
Purification notes	Purified from ascites by protein A affinity chromatography.
Clonality	Monoclonal
Clone number	IVA30
Isotype	IgG2

Applications

Our [Abpromise guarantee](#) covers the use of **ab3919** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent dilution.
IP		Use at an assay dependent dilution.

Application notes Is unsuitable for WB.

Target

Cellular localization Cell Membrane

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors