

Product datasheet

Anti-SOX15 antibody ab55992

1 Image

Overview

Product name	Anti-SOX15 antibody
Description	Chicken polyclonal to SOX15
Host species	Chicken
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment, corresponding to amino acids 130-233 of Human SOX15

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
Storage buffer	Preservative: None Constituents: PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgY

Applications

Our [Abpromise guarantee](#) covers the use of **ab55992** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		1/2000. Predicted molecular weight: 25 kDa.

Target

Function	Binds to the 5'-AACAAAT-3' sequence.
-----------------	--------------------------------------

Tissue specificity	Widely expressed in fetal and adult tissues examined, highest level found in fetal spinal cord and adult brain and testis.
Sequence similarities	Contains 1 HMG box DNA-binding domain.
Cellular localization	Nucleus.

Images



Anti-SOX15 antibody (ab55992) at 1/2000 dilution + E coli-derived fusion protein as test antigen.

Secondary

Goat anti-IgY-HRP dilution at 1/1000 dilution

Predicted band size: 25 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors