

Product datasheet

Anti-DZIP3 antibody ab56657

★★★★☆ 2 Abreviews 1 References 1 Image

Overview

Product name	Anti-DZIP3 antibody
Description	Mouse monoclonal to DZIP3
Host species	Mouse
Tested applications	Suitable for: IP, WB
Species reactivity	Reacts with: Recombinant fragment Predicted to work with: Human
Immunogen	Recombinant fragment: PSAGLRSDPS IMNWERITDR LKTAFPQQTR KELTDFLRKL KDAYGKSLSE LTFDEMCKI SQFIDPKKSQ SQGKSVSNVN CVSPSHSPSQ PDAAQPPKPA , corresponding to amino acids 1021-1121 of Human DZIP3 Run BLAST with ExPASy Run BLAST with NCBI

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None PBS, pH 7.2
Purity	Protein G purified
Clonality	Monoclonal
Isotype	IgG2a
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab56657** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IP	★★★★☆	Use at an assay dependent concentration.


Application	Abreviews	Notes
WB		Use a concentration of 1 - 5 µg/ml. This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

Target

Relevance E3 Ubiquitin ligase proteins mediate ubiquitination and subsequent proteasomal degradation of target proteins. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Able to specifically bind RNA.

Cellular localization Cytoplasm

Images



Western blot against tagged recombinant protein immunogen using ab56657 DZIP3 antibody at 1 µg/ml. Predicted band size of immunogen is 37 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors