

Product datasheet

Anti-rSec6 antibody ab56979

[1 References](#) [1 Image](#)

Overview

Product name	Anti-rSec6 antibody
Description	Mouse monoclonal to rSec6
Host species	Mouse
Tested applications	Suitable for: WB, ICC/IF
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment: SGFGEDVDGY CDTIVAVAEV IKLTDPSELLY LEVSTLVSKY PDIRDDHIGA LLAVRGDASR DMKQTIMETL EQGPAQASPS YVPLFKDIVV PSLNVAKLLK , corresponding to amino acids 646-746 of Human rSec6 Run BLAST with ExPASy Run BLAST with NCBI

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None PBS, pH 7.2
Purity	Protein G purified
Clonality	Monoclonal
Isotype	IgG2b
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab56979** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 - 5 µg/ml.


Application	Abreviews	Notes
ICC/IF		Use at an assay dependent concentration. PubMed: 24109596

Target

Function Component of the exocyst complex involved in the docking of exocytic vesicles with fusion sites on the plasma membrane.

Sequence similarities Belongs to the SEC6 family.

Images



Western blot against tagged recombinant protein immunogen using ab56979 rSec6 antibody at 1ug/ml. Predicted band size of immunogen is 37 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors