

Product datasheet

Anti-PCTP antibody ab57479

★★★★☆ 3 Abreviews 1 Image

Overview

Product name	Anti-PCTP antibody
Description	Mouse monoclonal to PCTP
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Mouse, Human
Immunogen	Recombinant full length protein, corresponding to amino acids 1-215 of Human PCTP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None PBS, pH 7.2
Purity	Protein G purified
Clonality	Monoclonal
Isotype	IgG2b
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab57479** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB	★★★★☆	Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 25 kDa.

Target

Function	Catalyzes the transfer of phosphatidylcholine between membranes. Binds a single lipid molecule.
Tissue specificity	Highest expression in liver, placenta, testis, kidney and heart. Low levels in brain and lung. No expression detected in thymus.
Sequence similarities	Contains 1 START domain.
Cellular localization	Cytoplasm.

Images



PCTP antibody (ab57479) at 1ug/lane + K-562 cell lysate at 25ug/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors