

Product datasheet

Mouse antibody [B-Z1] - Azide free ab59495

Overview

Product name	Mouse antibody [B-Z1] - Azide free
Tested applications	Suitable for: Flow Cyt
Immunogen	Rat monocytes and macrophages

Properties

Form	Liquid
Storage instructions	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: PBS
Purity	Ion Exchange Chromatography
Purification notes	ab59495 is sterile-filtered through 0.22 µm and treated to remove endotoxins.
Clonality	Monoclonal
Clone number	B-Z1
Myeloma	x63-Ag8.653
Isotype	IgG1

Applications

Our [Abpromise guarantee](#) covers the use of **ab59495** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		

Application notes	Flow Cyt: Use at an assay dependent dilution. Negative control for biological test. Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.
--------------------------	--

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors