

## Product datasheet

# Anti-MFAP3L antibody ab62188

1 Image

### Overview

<b>Product name</b>	Anti-MFAP3L antibody
<b>Description</b>	Rabbit polyclonal to MFAP3L
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>Species reactivity</b>	<b>Reacts with:</b> Human <b>Predicted to work with:</b> Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog
<b>Immunogen</b>	A 14-22 amino acid region within synthetic peptide: MHDSGLLNIT KVSFSDRGKY TCVASNIYGT VNNTVTLRVI FTSGDMGVYY, corresponding to N terminal amino acids 1-50 of Human MFAP3L <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
<b>Positive control</b>	Jurkat cell lysate.

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>Storage buffer</b>	Preservative: None Constituents: 2% Sucrose, PBS
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab62188** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 2.5 µg/ml. Detects a band of approximately 45 kDa (predicted molecular weight: 45 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## Target


---

<b>Tissue specificity</b>	Highly expressed in testis.
<b>Sequence similarities</b>	Contains 1 Ig-like C2-type (immunoglobulin-like) domain.
<b>Cellular localization</b>	Cell membrane.

---

## Images

---



**Lane 1 :** Marker

**Lane 2 :** Anti-MFAP3L antibody (ab62188) at 2.5 µg/ml

**Lane 2 :** Jurkat cell lysate at 10 µg

### Secondary

**Lane 2 :** HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 45 kDa

**Observed band size:** 45 kDa

Gel concentration 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors