

Product datasheet

Anti-CFHL2 antibody ab66976

1 Image

Overview

Product name	Anti-CFHL2 antibody
Description	Mouse polyclonal to CFHL2
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Full length human CFHL2 protein

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 1X PBS, pH 7.2
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab66976** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent dilution. Detects a band of approximately 31 kDa (predicted molecular weight: 31 kDa).

Target

Function Might be involved in complement regulation. Can associate with lipoproteins and may play a role

in lipid metabolism.

Tissue specificity

Expressed by the liver and secreted in plasma.

Sequence similarities

Contains 4 Sushi (CCP/SCR) domains.


Post-translational modifications

N-glycosylated.

Cellular localization

Secreted.

Images



All lanes : Anti-CFHL2 antibody (ab66976)

Lane 1 : CFHL2 transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Predicted band size: 31 kDa

Observed band size: 31 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors