

Product datasheet

Anti-CHIT1 antibody ab72574

[1 References](#) [1 Image](#)

Overview

Product name	Anti-CHIT1 antibody
Description	Mouse polyclonal to CHIT1
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Full length human protein (AAI03696).
Positive control	CHIT1 transfected 293T cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 1X PBS, pH 7.2
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab72574** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Detects a band of approximately 52 kDa (predicted molecular weight: 52 kDa).

Target

Function	Degrades chitin, chitotriose and chitobiose. May participate in the defense against nematodes and other pathogens. Isoform 3 has no enzymatic activity.
Tissue specificity	Detected in spleen. Secreted by cultured macrophages.
Sequence similarities	Belongs to the glycosyl hydrolase 18 family. Chitinase class II subfamily. Contains 1 chitin-binding type-2 domain.
Cellular localization	Secreted. Lysosome. A small proportion is lysosomal.

Images



All lanes : Anti-CHIT1 antibody (ab72574) at 1/500 dilution

Lane 1 : CHIT1 transfected 293T cell lysate

Lane 2 : non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse (H&L)-HRP Conjugated IgG at 1/2500 dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa

Blocking in 5% milk in PBST.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors