


## Product datasheet

### Anti-ZNF555 antibody ab75496

1 Image

#### Overview

Product name	Anti-ZNF555 antibody
Description	Rabbit polyclonal to ZNF555
Host species	Rabbit
Tested applications	<b>Suitable for:</b> WB
Species reactivity	<b>Reacts with:</b> Human <b>Predicted to work with:</b> Horse, Cow, Cat, Dog, Pig 
Immunogen	A synthetic peptide corresponding to a region of Human ZNF555
Positive control	Fetal Heart lysate.

#### Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

#### Applications


**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab75496 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 73 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## Target

<b>Function</b>	May be involved in transcriptional regulation.
<b>Sequence similarities</b>	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 15 C2H2-type zinc fingers. Contains 1 KRAB domain.
<b>Cellular localization</b>	Nucleus.

## Images



Anti-ZNF555 antibody (ab75496) at 1 µg/ml (in 5% skim milk / PBS buffer) + Fetal heart lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 73 kDa

**Observed band size:** 73 kDa

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors