

Product datasheet

Anti-KRTHA2 antibody ab76819

2 Images

Overview

Product name	Anti-KRTHA2 antibody
Description	Mouse polyclonal to KRTHA2
Host species	Mouse
Tested applications	Suitable for: ELISA, WB
Species reactivity	Reacts with: Human
Immunogen	Full length Human protein
Positive control	Lysate from Human placenta

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: PBS, pH 7.2
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab76819** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.
WB		1/500 - 1/1000. Predicted molecular weight: 50 kDa.

Target


Tissue specificity

Restricted to the hair cuticle.

Sequence similarities

Belongs to the intermediate filament family.

Images



Western blot - Anti-KRTHA2 antibody (ab76819)


Anti-KRTHA2 antibody (ab76819) at 1/500 dilution + Human placenta lysate at 25 μ g

Secondary

Goat Anti-Mouse IgG (H&L) HRP-conjugated at 1/2500 dilution

Predicted band size: 50 kDa

Observed band size: 48 kDa



Western blot - Anti-KRTHA2 antibody (ab76819)

All lanes : Anti-KRTHA2 antibody (ab76819) at 1/500 dilution

Lane 1 : KRTHA2 transfected 293T cell line

Lane 2 : non transfected 293T cell line

Lysates/proteins at 25 μ g per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L) HRP-conjugated at 1/2500 dilution

Predicted band size: 50 kDa

Observed band size: 50 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors