

Product datasheet

Anti-METTL7A antibody ab80235

1 Image

Overview

Product name	Anti-METTL7A antibody
Description	Rabbit polyclonal to METTL7A
Host species	Rabbit
Tested applications	Suitable for: IHC-P, ELISA
Species reactivity	Reacts with: Human
Immunogen	KLH conjugated synthetic peptide selected from the center region of human METTL7A.
Positive control	Human hepatocarcinoma tissue.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: 0.09% Sodium Azide Constituents: PBS
Purity	Protein A purified
Purification notes	ab80235 is purified through a protein A column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab80235** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/10 - 1/50.
ELISA		1/1000.

Target

Function	Probable methyltransferase.
Sequence similarities	Belongs to the methyltransferase superfamily.
Cellular localization	Lipid droplet. Endoplasmic reticulum. Membrane.

Images



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-METTL7A antibody (ab80235)

ab80235 at 1/10 dilution, staining METTL7A in human hepatocarcinoma by Immunohistochemistry using formalin-fixed, paraffin-embedded tissue, followed by a peroxidase-conjugated secondary antibody and DAB staining.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors