

Product datasheet

Anti-LYSMD4 antibody ab89035

2 Images

Overview

<b>Product name</b>	Anti-LYSMD4 antibody
<b>Description</b>	Mouse polyclonal to LYSMD4
<b>Host species</b>	Mouse
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	Full length protein, corresponding to amino acids 1-297 of Human LYSMD4 (NP_689662)
<b>Positive control</b>	Kidney lysate

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Storage buffer</b>	Constituent: Whole serum
<b>Purity</b>	Whole antiserum
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab89035** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Predicted molecular weight: 32 kDa.

Target


**Relevance** The specific function of LYSMD4 is not yet known.

## Images



Anti-LYSMD4 antibody (ab89035) at 1/500 dilution + Human kidney lysate at 50  $\mu$ g

**Predicted band size:** 32 kDa



**All lanes :** Anti-LYSMD4 antibody (ab89035) at 1/500 dilution

**Lane 1 :** LYSMD4-transfected 293T cell lysate

**Lane 2 :** non transfected lysate

Lysates/proteins at 25  $\mu$ g per lane.

**Predicted band size:** 32 kDa

**Observed band size:** 32 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors