

Product datasheet

Anti-DNAJC10 antibody ab98044

1 Image

Overview

Product name	Anti-DNAJC10 antibody
Description	Rabbit polyclonal to DNAJC10
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 35-84 (DFYSLGVSQ TASSREIRQA FKKLALKLHP DKNPNNPNAH GDFLKINRAY) of Human DNAJC10 (NP_061854). Run BLAST with ExPASy Run BLAST with NCBI
Positive control	ACHN cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 2% Sucrose, PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab98044** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 91 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

Target

Relevance	DNAJC10 (DnaJ homolog, subfamily C, member 10) is a co-chaperone which is thought to have a role in protein folding and translocation across the endoplasmic reticulum membrane. It may act as a co-chaperone for HSPA5.
Cellular localization	Secreted. Endoplasmic reticulum lumen.

Images



Anti-DNAJC10 antibody (ab98044) at 1 µg/ml

+ ACHN cell lysate at 10 µg

Predicted band size: 91 kDa

Gel concentration: 6-18%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors