

## Product datasheet

# Mouse monoclonal [L-1] Anti-Chicken lambda light chain ab99730

### Overview

<b>Product name</b>	Mouse monoclonal [L-1] Anti-Chicken lambda light chain
<b>Host species</b>	Mouse
<b>Target species</b>	Chicken
<b>Tested applications</b>	<b>Suitable for:</b> ICC/IF, IHC-Fr, ELISA, SDS-PAGE, Flow Cyt
<b>Immunogen</b>	The details of the immunogen for this antibody are not available.

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 8.20 Constituent: 100% Borate buffered saline  No preservatives or amine-containing buffer salts added
<b>Purity</b>	Protein G purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	L-1
<b>Isotype</b>	IgG1
<b>Light chain type</b>	kappa

### Applications

Our [Abpromise guarantee](#) covers the use of **ab99730** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration. Acetone fixed sections

Application	Abreviews	Notes
ELISA		Use at an assay dependent concentration.
SDS-PAGE		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. <a href="#">ab170190</a> - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors