

Product datasheet

Rat monoclonal [SB76b]Anti-Mouse IgG3 H&L ab99838

Overview

Product name	Rat monoclonal [SB76b]Anti-Mouse IgG3 H&L
Description	Rat monoclonal [SB76b]Anti-Mouse IgG3 H&L
Host species	Rat
Target species	Mouse
Specificity	Reacts with Mouse IgG3 as demonstrated by ELISA. Minimal reaction with other Mouse isotypes. May react with IgG from other species.
Tested applications	Suitable for: Flow Cyt, ELISA

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 8.20 Constituent: 100% Borate buffered saline
	No preservatives or amine-containing buffer salts added.
Purity	Protein A purified
Clonality	Monoclonal
Clone number	SB76b
Isotype	IgG2a

Applications

Our [Abpromise guarantee](#) covers the use of **ab99838** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use a concentration of 1 - 10 µg/ml.
ELISA		Use 1µg for 10 ⁶ cells.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors