

Product datasheet

Anti-ANKLE1 antibody ab222029

1 Image

Overview

Product name	Anti-ANKLE1 antibody
Description	Rabbit polyclonal to ANKLE1
Host species	Rabbit

Tested Applications & Species

Application	Species
IHC-P	Human

See all applications and species data

Immunogen Recombinant fragment corresponding to Human ANKLE1 aa 350-450.
Database link: [Q8NAG6](#)

[Run BLAST with](#)

[Run BLAST with](#)

Positive control IHC-P: Human bone marrow.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of **ab222029** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors