

Product datasheet

APC Mouse IgG2b [PLPV219] - Isotype Control ab91534

Overview

Product name	APC Mouse IgG2b [PLPV219] - Isotype Control
Conjugation	APC. Ex: 645nm, Em: 660nm
Specificity	ab91534 was produced against a synthetic hapten, which is normally not present in Humans or animals. Therefore, ab91534 does not react with any antigen of Human or animal cells.
Tested applications	Suitable for: Flow Cyt
Immunogen	Vegetal protein
General notes	25µg is sufficient for 50 tests

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Store In the Dark.
Storage buffer	pH: 7.20 Preservative: 0.1% Sodium azide Constituents: PBS, 1% BSA
Purity	Protein G purified
Clonality	Monoclonal
Clone number	PLPV219
Isotype	IgG2b

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of **ab91534** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		

Application notes Flow Cyt: Use 20µl for 10⁶ cells.
ab91534 is a negative isotype control

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors