

Product datasheet

Beagle Complement Preserved Serum ab155108

Overview

Product name	Beagle Complement Preserved Serum
Host species	Dog
General notes	Beagle Complement serum is collected from spontaneously clotted whole blood and centrifuged to separate the serum from the clot. Beagle Complement is ideal for hemolytic assay where hemolysis of IgG antibodies are required. Complement is preservative free. Raised in: Dog Purity: Whole antiserum

Properties

Form	Liquid
Storage instructions	Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle.
Purity	Whole antiserum

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors