

Product datasheet

Biotin Anti-C3 α /C3 α (desArg)/C3 antibody [D17/1] ab121050

[1 Image](#)

Overview

Product name	Biotin Anti-C3 α /C3 α (desArg)/C3 antibody [D17/1]
Description	Biotin Mouse monoclonal [D17/1] to C3 α /C3 α (desArg)/C3
Host species	Mouse
Conjugation	Biotin
Tested applications	Suitable for: IHC-P
Species reactivity	Reacts with: Human
Immunogen	Human C3 α /C3 α des Arg.
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties


Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.40 Preservative: 0.1% Sodium azide Constituents: 0.82% Sodium chloride, 99% Phosphate Buffer
Purity	Protein G purified
Purification notes	Protein-A/G purified.
Clonality	Monoclonal
Clone number	D17/1
Isotype	IgG1

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab121050 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent concentration.

Images



Immunohistochemistry of human appendix section staining C3a/C3a(desArg)/C3 with ab121050 (1:30)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Biotin Anti-C3a/C3a(desArg)/C3 antibody [D17/1] (ab121050)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors