

Product datasheet

Biotin Anti-Integrin alpha 5 antibody [5H10-27]  
ab25189

1 Image

Overview

<b>Product name</b>	Biotin Anti-Integrin alpha 5 antibody [5H10-27]
<b>Description</b>	Biotin Rat monoclonal [5H10-27] to Integrin alpha 5
<b>Host species</b>	Rat
<b>Conjugation</b>	Biotin
<b>Specificity</b>	ab25189 recognises alpha 5 subunit of the integrin alpha 5 beta 1 fibronectin receptor.
<b>Tested applications</b>	<b>Suitable for:</b> Flow Cyt
<b>Immunogen</b>	Tissue, cells or virus corresponding to Mouse Integrin alpha 5. C57BL/6 x A/J F1 mouse mast cell line MC/9

General notes

The monoclonal antibody 5H10-27 reacts with the alpha 5 chain (VLA-5) of the CD49e/CD29 heterodimeric fibronectin receptor. Soluble 5H10-27 antibody inhibits VLA-5-mediated functions in vitro. Immobilized 5H10-27 has also been demonstrated to co-stimulate the proliferative response of CD8+ T cells to plate-bound anti-CD3 monoclonal antibody.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS

<b>Purity</b>	Affinity purified
<b>Primary antibody notes</b>	The monoclonal antibody 5H10-27 reacts with the alpha 5 chain (VLA-5) of the CD49e/CD29 heterodimeric fibronectin receptor. Soluble 5H10-27 antibody inhibits VLA-5-mediated functions in vitro. Immobilized 5H10-27 has also been demonstrated to co-stimulate the proliferative response of CD8+ T cells to plate-bound anti-CD3 monoclonal antibody.
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	5H10-27
<b>Isotype</b>	IgG2a
<b>Light chain type</b>	kappa

## Applications


**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab25189 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use 1µg for 10 <sup>6</sup> cells. <a href="#">ab18445</a> - Rat monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

## Target

<b>Function</b>	Integrin alpha-5/beta-1 is a receptor for fibronectin and fibrinogen. It recognizes the sequence R-G-D in its ligands. In case of HIV-1 infection, the interaction with extracellular viral Tat protein seems to enhance angiogenesis in Kaposi's sarcoma lesions.
<b>Sequence similarities</b>	Belongs to the integrin alpha chain family. Contains 7 FG-GAP repeats.
<b>Cellular localization</b>	Membrane.

## Images



Flow Cytometry - Biotin Anti-Integrin alpha 5 antibody [5H10-27] (ab25189)

Flow cytometry analysis staining Integrin alpha 5 in cells using ab25189 at a dilution of 1 µg/10<sup>6</sup> cells.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors