

Product datasheet

Bovine Heart Mitochondria ab110338

3 References

Overview

Product name Bovine Heart Mitochondria

General notes

2 mg of bovine heart mitochondria in iso-osmotic buffer. This sample can act as a positive control in immunocapture applications and Blue Native PAGE.

1 mg solubilized mitochondria for Immunocapture experiments. For the Immunocapture of OXPHOS complexes, it is recommended that 180 µl of the mitochondrial membranes supplied be solubilized by the addition of 20 µl of 10% lauryl maltoside. Insoluble material is removed by centrifugation at 72,000g before immunocapture of complexes from the supernatant. For more detail see Abcam's suggested protocols.

Tested applications **Suitable for:** IP, SDS-PAGE, WB

Properties

Mycoplasma free Yes

Form Liquid

Storage instructions Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

Storage buffer Constituents: 0.02% PMSF, 0.16% Tris, 0.89% EDTA, 8.56% Sucrose

Lysate notes 2 mg of bovine heart mitochondria in iso-osmotic buffer. This sample can act as a positive control in immunocapture applications and Blue Native PAGE.

1 mg solubilized mitochondria for Immunocapture experiments. For the Immunocapture of OXPHOS complexes, it is recommended that 180 µl of the mitochondrial membranes supplied be solubilized by the addition of 20 µl of 10% lauryl maltoside. Insoluble material is removed by centrifugation at 72,000g before immunocapture of complexes from the supernatant. For more detail see Abcam's suggested protocols.

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab110338 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IP		Use at an assay dependent concentration.
SDS-PAGE		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors