



Product datasheet

Anti-BROX antibody - C-terminal ab193008

1 Image

Overview

Product name	Anti-BROX antibody - C-terminal
Description	Rabbit polyclonal to BROX - C-terminal
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	<p>Synthetic peptide within Human BROX aa 360-409 (C terminal). The exact sequence is proprietary. NCBI Accession No. NP_653296.</p> <p>Sequence:</p> <p>PTSVQWTPETLAAFDLTKRPKDDSTKPKPEEEVKPKV EPDIKPQKDTGCV</p> <p>Database link: Q5VW32</p> <p style="text-align: right;">  Run BLAST with  Run BLAST with </p>
Positive control	Jurkat whole cell lysate.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	<p>Preservative: 0.09% Sodium azide</p> <p>Constituents: 2% Sucrose, PBS</p>
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab193008** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 46 kDa.

Target

Sequence similarities	Belongs to the BROX family. Contains 1 BRO1 domain.
Cellular localization	Membrane.

Images



Anti-BROX antibody - C-terminal (ab193008) at 1 µg/ml + Jurkat whole cell lysate at 10 µg

Predicted band size: 46 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors