

Product datasheet

Anti-C13orf30 antibody ab243593

1 Image

Overview

Product name	Anti-C13orf30 antibody
Description	Rabbit polyclonal to C13orf30
Host species	Rabbit

Tested Applications & Species

Application	Species
IHC-P	Human

See all applications and species data

Immunogen Recombinant fragment corresponding to Human C13orf30 aa 1 to the C-terminus.
Database link: [Q8N7L0](#)

[Run BLAST with](#)

[Run BLAST with](#)

Positive control IHC-P: Human colon tissue.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of **ab243593** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

✔ Guaranteed

Tested applications are guaranteed to work and covered by our Abpromise guarantee.

⚠ Predicted

Predicted to work for this combination of applications and species but not guaranteed.

✘ Incompatible

Does not work for this combination of applications and species.


Application	Species
IHC-P	✔ Human

Application	Abreviews	Notes
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

Target

Sequence similarities Belongs to the FAM216 family.

Images



Paraffin-embedded human colon tissue stained for C13orf30 using ab243593 at 1/50 dilution in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C13orf30 antibody (ab243593)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors