


Product datasheet

Anti-C7ORF50 antibody ab214080

1 Image

Overview

Product name	Anti-C7ORF50 antibody
Description	Rabbit polyclonal to C7ORF50
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human Predicted to work with: Rat 
Immunogen	Synthetic peptide corresponding to Human C7ORF50 aa 151-165. Sequence: LQGRARELTVQKAEA Database link: Q9BRJ6 Run BLAST with Run BLAST with
Positive control	Lysate of HEK293T cells overexpressing C7ORF50.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.40 Preservative: 0.098% Sodium azide Constituent: 99% PBS
Purity	IgG fraction
Clonality	Polyclonal
Isotype	IgG

Applications


Our [Abpromise guarantee](#) covers the use of **ab214080** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/2000. Predicted molecular weight: 22 kDa.

Target

Images



Anti-C7ORF50 antibody (ab214080) at 1/1000 dilution + Lysate of HEK293T cells overexpressing C7ORF50

Secondary

goat anti-rabbit IgG peroxidase-conjugated

Developed using the ECL technique.

Predicted band size: 22 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors