

Product datasheet

Anti-CD134 / OX40L receptor antibody [Ber-ACT35]
ab234239

1 Image

Overview

Product name	Anti-CD134 / OX40L receptor antibody [Ber-ACT35]
Description	Mouse monoclonal [Ber-ACT35] to CD134 / OX40L receptor
Host species	Mouse
Tested applications	Suitable for: Flow Cyt
Species reactivity	Reacts with: Human Predicted to work with: Non human primates ▲
Immunogen	Tissue, cells or virus corresponding to Human CD134/ OX40L receptor. (HTLV 1-transformed HUT-102 cells).
Positive control	Flow Cyt: Human PHA-activated leukocytes.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.40 Preservative: 0.0975% Sodium azide Constituent: PBS
Purity	Protein A purified
Purification notes	Purity > 95% (by SDS-PAGE).
Clonality	Monoclonal
Clone number	Ber-ACT35
Isotype	IgG1

Applications

Our [Abpromise guarantee](#) covers the use of **ab234239** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent concentration.

Target

Function	Receptor for TNFSF4/OX40L/GP34.
Sequence similarities	Contains 4 TNFR-Cys repeats.
Cellular localization	Membrane.

Images



Flow cytometric analysis of human PHA-activated leukocytes labeling CD134 / OX40L receptor with ab234239, followed by Goat anti-mouse APC (red) compared with negative controls. Surface staining.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors