

Product datasheet

Anti-CD147 antibody [EPR18008-8] ab188190

Recombinant RabMAb

[2 References](#) [9 Images](#)

Overview

Product name	Anti-CD147 antibody [EPR18008-8]
Description	Rabbit monoclonal [EPR18008-8] to CD147
Host species	Rabbit
Tested applications	Suitable for: WB, IHC-P, ICC/IF, Flow Cyt, IP, Sandwich ELISA
Species reactivity	Reacts with: Mouse
Immunogen	Recombinant fragment within Mouse CD147 aa 100-350. The exact sequence is proprietary. Database link: P18572
Positive control	WB: Mouse CD147 recombinant protein fragment; RAW 264.7, WEHI-3 and bEND.3 whole cell lysates; mouse liver lysate. IHC-P: Mouse intestine tissue. ICC/IF: WEHI-231 and bEND.3 cells. Flow Cyt: Mouse thymocytes. IP: RAW 264.7 whole cell lysate.
General notes	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p> <p>Reproducibility is key to advancing scientific discovery and accelerating scientists' next breakthrough.</p> <p>Abcam is leading the way with our range of recombinant antibodies, knockout-validated antibodies and knockout cell lines, all of which support improved reproducibility.</p> <p>We are also planning to innovate the way in which we present recommended applications and species on our product datasheets, so that only applications & species that have been tested in our own labs, our suppliers or by selected trusted collaborators are covered by our Abpromise[™] guarantee.</p> <p>In preparation for this, we have started to update the applications & species that this product is Abpromise guaranteed for.</p> <p>We are also updating the applications & species that this product has been "predicted to work with," however this information is not covered by our Abpromise guarantee.</p>

Applications & species from publications and Abreviews that have not been tested in our own labs or in those of our suppliers are not covered by the Abpromise guarantee.

Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, as well as customer reviews and Q&As.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol, 0.05% BSA
Purity	Protein A purified
Clonality	Monoclonal
Clone number	EPR18008-8
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab188190** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/5000. Detects a band of approximately 55 kDa (predicted molecular weight: 42 kDa).
IHC-P		1/8000. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.
ICC/IF		1/250.
Flow Cyt		1/200. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
IP		1/50.
Sandwich ELISA		Use at an assay dependent concentration.

Target

Function	Plays pivotal roles in spermatogenesis, embryo implantation, neural network formation and tumor
-----------------	---

progression. Stimulates adjacent fibroblasts to produce matrix metalloproteinases (MMPs). May target monocarboxylate transporters SLC16A1, SLC16A3 and SLC16A8 to plasma membranes of retinal pigment epithelium and neural retina. Seems to be a receptor for oligomannosidic glycans. In vitro, promotes outgrowth of astrocytic processes.

Tissue specificity

Present only in vascular endothelium in non-neoplastic regions of the brain, whereas it is present in tumor cells but not in proliferating blood vessels in malignant gliomas.

Sequence similarities

Contains 1 Ig-like C2-type (immunoglobulin-like) domain.
Contains 1 Ig-like V-type (immunoglobulin-like) domain.


Post-translational modifications

N-glycosylated.

Cellular localization

Cell membrane. Melanosome. Colocalizes with SLC16A1 and SLC16A8 (By similarity). Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

Images



All lanes : Anti-CD147 antibody [EPR18008-8] (ab188190) at 1/5000 dilution

Lane 1 : RAW 264.7 (Mouse macrophage cell line transformed with Abelson murine leukemia virus) whole cell lysate

Lane 2 : Mouse liver lysate

Lane 3 : WEHI-3 (Mouse leukemia cell line) whole cell lysate

Lane 4 : bEnd.3 (Mouse brain endothelioma cell line) whole cell lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/100000 dilution

Predicted band size: 42 kDa


Observed band size: 45-55 kDa

[why is the actual band size different from the predicted?](#)

Blocking/Dilution buffer: 5% NFDm/TBST.

Exposure times: Lane 1/2: 3 minutes; Lane 3: 10 seconds; Lane 4: 1 second.

The expression profile observed is consistent with what has been described in the literature (PMID: 16721788; 23966157).



Immunocytochemistry/ Immunofluorescence - Anti-CD147 antibody [EPR18008-8] (ab188190)


Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized bEND.3 (Mouse brain capillary endothelial cell line) cells labeling CD147 with ab188190 at 1/250 dilution, followed by by Goat Anti-Rabbit IgG H&L (Alexa Fluor[®] 488) (ab150077) secondary antibody at 1/1000 dilution (green). Confocal image showing cell membrane staining on bEND.3 cell line. The nuclear counter stain is DAPI (blue).

Tubulin is detected with Anti-alpha Tubulin mouse MAb (ab7291) at 1/1000 dilution and Goat Anti-Mouse IgG H&L (Alexa Fluor[®] 594) (ab150120) secondary antibody at 1/1000 dilution (red).

The negative controls are as follows:

-ve control 1: ab188190 at 1/250 dilution followed by Goat Anti-Mouse IgG H&L (Alexa Fluor[®] 594) (ab150120) secondary antibody at 1/1000 dilution.

-ve control 2: Anti-alpha Tubulin mouse MAb (ab7291) at 1/1000 dilution followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor[®] 488) (ab150077) secondary antibody at 1/1000 dilution.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD147 antibody [EPR18008-8] (ab188190)

Immunohistochemical analysis of paraffin-embedded mouse intestine tissue labeling CD147 with ab188190 at 1/8000 dilution, followed by Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/500 dilution. Membrane staining on mouse intestine is observed. Counter stained with Hematoxylin.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/500 dilution.

Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.



Immunoprecipitation - Anti-CD147 antibody [EPR18008-8] (ab188190)

CD147 was immunoprecipitated from 1 mg of RAW 264.7 (Mouse macrophage cell line transformed with Abelson murine leukemia virus) whole cell lysate with ab188190 at 1/50 dilution. Western blot was performed from the immunoprecipitate using ab188190 at 1/5000 dilution. VeriBlot for IP Detection Reagent (HRP) (ab131366), was used for detection at 1/10000 dilution.


Lane 1: RAW 264.7 whole cell lysate 10µg (Input).

Lane 2: ab188190 IP in RAW 264.7 whole cell lysate.

Lane 3: Rabbit monoclonal IgG (ab172730) instead of ab188190 in RAW 264.7 whole cell lysate.


Blocking and dilution buffer and concentration: 5% NFDm/TBST.

Exposure time: 1 second.



Flow Cytometry - Anti-CD147 antibody [EPR18008-8] (ab188190)

Flow cytometric analysis of fresh mouse thymocytes labeling CD147 with ab188190 at 1/200 dilution (red) compared with Rabbit IgG, monoclonal [EPR25A] - Isotype Control (ab172730; black) and an unlabeled control (cells without incubation with primary antibody and secondary antibody; blue). Goat Anti-Rabbit IgG (Alexa Fluor® 488) at 1/500 dilution was used as the secondary antibody.




Sandwich ELISA - Anti-CD147 antibody [EPR18008-8] (ab188190)

Standard Curve for CD147 (Analyte: Recombinant mouse CD147 protein) dilution range 0-500 pg/mL using Capture antibody at 0.2 ug/mL and Detector Antibody at 0.5 ug/mL. Secondary antibody: Peroxidase Streptavidin SA-HRP at 1/20000 dilution. Concentration of ab188190 may vary from lot to lot; please use this curve as guideline.

Washing buffer: 1X PBST

Blocking/Diluting buffer and concentration: 1% BSA/PBS



Immunocytochemistry/ Immunofluorescence - Anti-CD147 antibody [EPR18008-8] (ab188190)


Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized WEHI-231 (Mouse B Cell Lymphoma cell line) cells labeling CD147 with ab188190 at 1/250 dilution, followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor[®] 488) (ab150077) secondary antibody at 1/1000 dilution (green). Confocal image showing cell membrane staining on WEHI-231 cell line. The nuclear counter stain is DAPI (blue).

Tubulin is detected with Anti-alpha Tubulin mouse MAb (ab7291) at 1/1000 dilution and Goat Anti-Mouse IgG H&L (Alexa Fluor[®] 594) (ab150120) secondary antibody at 1/1000 dilution (red).

The negative controls are as follows:

-ve control 1: ab188190 at 1/250 dilution followed by Goat Anti-Mouse IgG H&L (Alexa Fluor[®] 594) (ab150120) secondary antibody at 1/1000 dilution.

-ve control 2: Anti-alpha Tubulin mouse MAb (ab7291) at 1/1000 dilution followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor[®] 488) (ab150077) secondary antibody at 1/1000 dilution.



Western blot - Anti-CD147 antibody [EPR18008-8] (ab188190)

Anti-CD147 antibody [EPR18008-8] (ab188190) at 1/5000 dilution + Mouse CD147 recombinant protein fragment at 0.002 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/100000 dilution

Predicted band size: 42 kDa


Observed band size: 28 kDa [why is the actual band size different from the predicted?](#)

Exposure time: 3 minutes

Blocking/Dilution buffer: 5% NFDM/TBST.

Recombinant protein fragment Mouse CD147 contains aa140-325 with a His-Tag®. It was made in house

Why choose a recombinant antibody?



- Research with confidence**
Consistent and reproducible results
- Long-term and scalable supply**
Recombinant technology
- Success from the first experiment**
Confirmed specificity
- Ethical standards compliant**
Animal-free production

Anti-CD147 antibody [EPR18008-8] (ab188190)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors