


Product datasheet

Anti-CD4 antibody [SK3] (PE/Cy7[®]) ab210372

1 Image

Overview

Product name	Anti-CD4 antibody [SK3] (PE/Cy7 [®])
Description	Mouse monoclonal [SK3] to CD4 (PE/Cy7 [®])
Host species	Mouse
Conjugation	PE/Cy7 [®] . Ex: 496nm, Em: 774nm
Tested applications	Suitable for: Flow Cyt
Species reactivity	Reacts with: Human Predicted to work with: Cynomolgus monkey, Rhesus monkey 
Immunogen	The details of the immunogen for this antibody are not available.
Positive control	Human peripheral blood lymphocytes
General notes	This antibody does not block binding of alternative clone RPA-T4, suggesting that they recognize different epitopes. This product or portions thereof is manufactured under license from Carnegie Mellon University under U.S. Patent Number 5, 268, 486 and related patents. Cy and CyDye are trademarks of GE Healthcare Limited.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Store In the Dark.
Storage buffer	pH: 7.2 Preservative: 0.09% Sodium azide Constituents: 0.12% Monobasic dihydrogen sodium phosphate, 0.87% Sodium chloride, 0.1% Gelatin
Purity	Affinity purified
Purification notes	ab210372 was purified from tissue culture supernatant via affinity chromatography. The purified antibody was conjugated under optimal conditions, with unreacted dye removed from the preparation.
Clonality	Monoclonal
Clone number	SK3
Isotype	IgG1

Light chain type

kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab210372** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use 5µl for 10 ⁵⁻⁸ cells. The antibody has been diluted for use at 5 µl (0.06 µg) per test, defined as the amount of antibody that will stain a cell sample in a final volume of approximately 100 µl.

Target

Function

Accessory protein for MHC class-II antigen/T-cell receptor interaction. May regulate T-cell activation. Induces the aggregation of lipid rafts.

Sequence similarities

Contains 3 Ig-like C2-type (immunoglobulin-like) domains.
Contains 1 Ig-like V-type (immunoglobulin-like) domain.


Post-translational modifications

Palmitoylation and association with LCK contribute to the enrichment of CD4 in lipid rafts.

Cellular localization

Cell membrane. Localizes to lipid rafts. Removed from plasma membrane by HIV-1 Nef protein that increases clathrin-dependent endocytosis of this antigen to target it to lysosomal degradation. Cell surface expression is also down-modulated by HIV-1 Envelope polyprotein gp160 that interacts with, and sequesters CD4 in the endoplasmic reticulum.

Images



Flow cytometric analysis of Human peripheral blood lymphocytes labeling CD4 with 5 µl (0.06 µg) ab210372 (solid line), compared with 0.06 µg PE/Cy7 Mouse IgG1 isotype control (dashed line).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors