

## Product datasheet

# Anti-Cdc6 antibody [Cdc6 9H8/5] ab20150

### 6 References

#### Overview

---

<b>Product name</b>	Anti-Cdc6 antibody [Cdc6 9H8/5]
<b>Description</b>	Mouse monoclonal [Cdc6 9H8/5] to Cdc6
<b>Host species</b>	Mouse
<b>Species reactivity</b>	<b>Reacts with:</b> <i>Saccharomyces cerevisiae</i>
<b>Immunogen</b>	Full length protein corresponding to <i>Saccharomyces cerevisiae</i> Cdc6.
<b>General notes</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Constituent: PBS
<b>Purity</b>	Protein A/G purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	Cdc6 9H8/5
<b>Myeloma</b>	unknown
<b>Isotype</b>	IgG1
<b>Light chain type</b>	unknown

#### Target

---

<b>Function</b>	Involved in the initiation of DNA replication. Also participates in checkpoint controls that ensure DNA replication is completed before mitosis is initiated.
-----------------	---

<b>Sequence similarities</b>	Belongs to the CDC6/cdc18 family.
<b>Cellular localization</b>	Nucleus. Cytoplasm. The protein is nuclear in G1 and cytoplasmic in S-phase cells.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors