

## Product datasheet

# Cell Lysis Buffer ab152163

21 References

### Overview

<b>Product name</b>	Cell Lysis Buffer
<b>Tested applications</b>	<b>Suitable for:</b> WB, IP
<b>General notes</b>	<p>ab152163 is a 10X Cell Lysis Buffer used to lyse cells under non-denaturing conditions.</p> <p>Recommended for preparing samples for use in western blotting and immunoprecipitation. For use with adherent or suspension cells.</p> <p>Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances.</p> <p>It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.</p>

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Store at -20°C. Stable for 12 months at -20°C
<b>Storage buffer</b>	<p>pH: 7.50</p> <p>Constituents: 0.29% Sodium EDTA, 1.12% Sodium pyrophosphate decahydrate, 8.8% Sodium chloride, 3.15% Tris HCl, 10% Triton-X-100, 0.38% EGTA, 0.001% Leupeptin, 0.19% Sodium orthovanadate, 0.216% Beta glycerophosphate</p>

### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab152163 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors