

Product datasheet

Anti-Chemerin antibody [Tigy-1] ab112520

★★★★☆ 1 Abreviews 1 References

Overview

Product name	Anti-Chemerin antibody [Tigy-1]
Description	Mouse monoclonal [Tigy-1] to Chemerin
Host species	Mouse
Tested applications	Suitable for: WB, ELISA
Species reactivity	Reacts with: Rat, Human
Immunogen	Recombinant fragment, corresponding to amino acids 23-157 of Rat Chemerin.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: 0.02% Sodium azide Constituent: 99% PBS
Purity	Protein G purified
Purification notes	ab112520 is purified from concentrated hybridoma tissue culture supernatant. Purity is above 95% (SDS-PAGE).
Clonality	Monoclonal
Clone number	Tigy-1
Isotype	IgG2a

Applications

Our [Abpromise guarantee](#) covers the use of **ab112520** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB	★★★★☆	1/1000. Predicted molecular weight: 19 kDa.
ELISA		Use at an assay dependent dilution.

Target

Tissue specificity	Highly expressed in skin (basal and suprabasal layers of the epidermis, hair follicles and endothelial cells). Also found in pancreas, liver, spleen, prostate, ovary, small intestine and colon.
Cellular localization	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors