

Product datasheet

Anti-Cysteine sulfonate antibody ab176487

1 Image

Overview

Product name	Anti-Cysteine sulfonate antibody
Description	Rabbit polyclonal to Cysteine sulfonate
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Species independent
Immunogen	Chemical/ Small Molecule corresponding to Cysteine sulfonate. Recombinant Human cysteine (sulfonate). Cysteamine backbone modified by dimedone with a five-carbon linker to KLH.
Positive control	Hela cell lysate and Hela heat shocked lysate.

Properties


Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.09% Sodium azide Constituents: 49% PBS, 50% Glycerol
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab176487** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000. ECL



All lanes : Anti-Cysteine sulfonate antibody (ab176487) at 1/1000 dilution

Lane 1 : HeLa cell lysate

Lane 2 : HeLa cell lysate, heat shocked

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors