

## Product datasheet

# Anti-DDDDK tag antibody - (Equivalent to FLAG antibodies from Sigma) ab122902

★☆☆☆☆ 1 Abreviews 1 Image

### Overview

<b>Product name</b>	Anti-DDDDK tag antibody - (Equivalent to FLAG antibodies from Sigma)
<b>Description</b>	Rabbit polyclonal to DDDDK tag - (Equivalent to FLAG antibodies from Sigma)
<b>Host species</b>	Rabbit
<b>Specificity</b>	<a href="#">ab12902</a> recognizes Flag-tagged fusion proteins.
<b>Tested applications</b>	<b>Suitable for:</b> WB, ELISA
<b>Species reactivity</b>	<b>Reacts with:</b> Species independent
<b>Immunogen</b>	Synthetic peptide corresponding to DDDDK tag conjugated to Keyhole Limpet Haemocyanin (KLH). Sequence: DYKDDDDK
	<a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
<b>Positive control</b>	Lysate of 293T cells transfected with the plasmid carrying the Flag-tagged PRMT6 gene.

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at -20°C.
<b>Storage buffer</b>	Preservative: 0.05% Sodium azide Constituent: 99% Whole serum
<b>Purity</b>	Whole antiserum
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab122902** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB	★☆☆☆☆	1/2000.
ELISA		Use at an assay dependent dilution.

## Target

**Relevance** This is a useful tool for the localisation and characterisation of DDDDK tagged proteins.

## Images



**All lanes** : Anti-DDDDK tag antibody - (Equivalent to FLAG antibodies from Sigma) (ab122902) at 1/2000 dilution

**Lane 1** : Lysate of 293T cells transfected with an empty vector

**Lane 2** : Lysate of 293T cells transfected with the plasmid carrying the Flag-tagged PRMT6 gene

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors