

Product datasheet

Donkey Anti-Goat IgG H&L (12nm Gold) preadsorbed ab105269

Overview

Product name	Donkey Anti-Goat IgG H&L (12nm Gold) preadsorbed	
Host species	Donkey	
Target species	Goat	
Tested applications	Suitable for: Electron Microscopy, WB	
Minimal cross-reactivity	Chicken, Guinea pig, Hamster, Horse, Human, Mouse, Rabbit, Rat	more details
Conjugation	Gold 12nm	

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 8.50 Preservative: 0.05% Sodium azide Constituents: 1.5% BSA, 0.87% Sodium chloride, 0.164% Sodium phosphate, 0.381% Sodium borate
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of **ab105269** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Electron Microscopy		1/20 - 1/40.
WB		1/20 - 1/40.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors