

Product datasheet

Donkey Anti-Goat IgG H&L (FITC) ab6881

[42 References](#) [2 Images](#)

Overview

Product name	Donkey Anti-Goat IgG H&L (FITC)
Host species	Donkey
Target species	Goat
Tested applications	Suitable for: IHC-P, IHC-Fr, ICC/IF, Immunomicroscopy, Flow Cyt, ELISA, WB
Immunogen	Full length native Goat IgG (purified).
Conjugation	FITC. Ex: 493nm, Em: 528nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.2 Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride, 1% BSA
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
Conjugation notes	Diaminotriazinylaminofluorescein (DTAF) (Molecular Weight 530 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 3.3 moles DTAF per mole of Rabbit IgG
Clonality	Polyclonal
Isotype	IgG


Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab6881 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ICC/IF		1/1000 - 1/5000.
Immunomicroscopy		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent dilution.
ELISA		1/10000 - 1/50000.
WB		1/5000 - 1/20000.

Images



Flow Cytometry - Donkey Anti-Goat IgG H&L (FITC)
(ab6881)

ab6881 was used at dilution 1/200 with the primary antibody [ab27985](#) in Flow Cyt. See the review on [ab27985](#).



Immunocytochemistry/ Immunofluorescence -
Donkey Anti-Goat IgG H&L (FITC) (ab6881)

ab6881 was used at dilution 1/500 with the primary antibody [ab26234](#) in ICC/IF. See the review on [ab26234](#).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors