

Product datasheet

Donkey Anti-Mouse IgG H&L (FITC) ab97029

★★★★★ [1 Abreviews](#) [5 References](#)

Overview

| | |
|----------------------------|---|
| Product name | Donkey Anti-Mouse IgG H&L (FITC) |
| Host species | Donkey |
| Target species | Mouse |
| Specificity | By immunoelectrophoresis and ELISA this antibody reacts specifically with mouse IgG and with light chains common to other mouse immunoglobulins. No antibody was detected against nonimmunoglobulin serum proteins. |
| Tested applications | Suitable for: ICC, IHC-P, Flow Cyt |
| Conjugation | FITC. Ex: 493nm, Em: 528nm |

Properties

| | |
|-----------------------------|--|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C. |
| Storage buffer | pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS |
| Purity | Immunogen affinity purified |
| Purification notes | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to FITC. |
| Conjugation notes | F/P ratio is 3-7 |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab97029 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|---------------|
| ICC | | 1/50 - 1/200. |
| IHC-P | | 1/50 - 1/500. |
| Flow Cyt | | 1/50 - 1/200. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors