

Product datasheet

Donkey anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed ab216778

1 References 1 Image

Overview

Product name	Donkey anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed
Description	Donkey anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed
Host species	Donkey
Target species	Mouse
Tested applications	Suitable for: WB, In-Cell ELISA
Minimal cross-reactivity	Chicken, Cow, Goat, Guinea pig, Horse, Human, Rabbit, Sheep, Syrian hamster more details
Immunogen	Full length protein corresponding to Mouse IgG. Whole molecule
Conjugation	IRDye® 680RD

Properties


Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Store In the Dark.
Storage buffer	Preservative: 0.01% Sodium azide Constituent: PBS Contains 10 mg/mL BSA (free of IgG and protease) as stabilizer.
Purity	Affinity purified
Purification notes	Affinity chromatography purified using antigens coupled to agarose.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of **ab216778** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/5000 - 1/25000.
In-Cell ELISA		1/200 - 1/800.

Images



Excitation and emission spectra of IRDye[®] 680RD Ex: 680, Em:694.

Spectrum - Donkey anti-Mouse IgG H&L (IRDye[®] 680RD) preadsorbed (ab216778)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors