

Product datasheet

Donkey Anti-Rabbit IgG H&L (Biotin) ab6801

12 References 2 Images

Overview

Product name	Donkey Anti-Rabbit IgG H&L (Biotin)
Host species	Donkey
Target species	Rabbit
Tested applications	Suitable for: WB, ELISA, Immunomicroscopy, Dot blot, ICC/IF, IHC-P, IHC-Fr
Immunogen	Rabbit IgG whole molecule
Conjugation	Biotin

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Potassium phosphate
Purity	Affinity purified
Purification notes	Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC). Biotin/Protein Ratio: 10-20 BAC molecules per Donkey IgG molecule.
Conjugation notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab6801 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		1/3000 - 1/15000.
ELISA		1/150000.
Immunomicroscopy		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
ICC/IF		1/500 - 1/3000.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

Images



Ab6801 was used at dilution 1/200 with the primary antibody [ab64869](#) in IHC-Fr. See the review on [ab64869](#).

Immunohistochemistry (Frozen sections) - Donkey
Anti-Rabbit IgG H&L (Biotin) (ab6801)



Ab6801 was used at dilution 1/200 with the primary antibody [ab39350](#) in IHC - Free Floating. See the review on [ab39350](#).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Donkey Anti-Rabbit IgG H&L (Biotin) (ab6801)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors