

Product datasheet

Donkey Anti-Sheep IgG H&L (Biotin) ab97123

2 References

Overview

Product name	Donkey Anti-Sheep IgG H&L (Biotin)
Description	Donkey polyclonal Secondary Antibody to Sheep IgG - H&L (Biotin)
Host species	Donkey
Target species	Sheep
Tested applications	Suitable for: ELISA, ICC/IF, IHC-P, WB
Immunogen	Other Immunogen Type corresponding to Sheep IgG.
Conjugation	Biotin

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
Purity	Immunogen affinity purified
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin. Biotinylated antibody was demonstrated with avidin/peroxidase when coated on microtiter wells.
Conjugation notes	Biotin/ antibody protein ratio is 4:1
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab97123** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/20000 - 1/400000. (primary)
ICC/IF		1/200 - 1/1000.
IHC-P		1/500 - 1/5000.
WB		1/20000 - 1/400000. (chemiluminescent)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors