

Product datasheet

Anti-EBNA 3C antibody ab16128

★★★★☆ 1 Abreviews 6 References

Overview

Product name	Anti-EBNA 3C antibody
Description	Sheep polyclonal to EBNA 3C
Host species	Sheep
Specificity	This antibody is specific for the Epstein Barr Virus Nuclear Antigen 3C.
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Species independent
Immunogen	Full length protein.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: 0.08% Sodium azide Constituent: PBS
Purity	Ammonium Sulphate Precipitation
Purification notes	Purified by ammonium sulfate precipitation This antibody is provided as a 0.2µm sterile filtered solution.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab16128** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB	★★★★☆	

Application notes WB: 1/200.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Target

Relevance

EBNA 3C is a latent viral nuclear protein expressed in EBV transformed lymphoblastic cell lines and is essential for EBV mediated transformation of these cell lines. This viral nuclear protein is essential for EBV mediated positive immunoblastic lymphomas, and functions as a transcription factor.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors