


Product datasheet

Anti-eIF5A2 + eIF5A antibody [EPR7412-50] - BSA and Azide free ab271915

Recombinant RabMAb

7 Images

Overview

Product name	Anti-eIF5A2 + eIF5A antibody [EPR7412-50] - BSA and Azide free
Description	Rabbit monoclonal [EPR7412-50] to eIF5A2 + eIF5A - BSA and Azide free
Host species	Rabbit
Tested applications	Suitable for: WB, IHC-P Unsuitable for: IP
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat 
Immunogen	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
Positive control	IHC-P: Recombinant Human eIF5A2 protein (ab99140) can be used as a positive control in WB. HepG2, HeLa, 293T and Raji cell lysates; Human colon adenocarcinoma and testis tissues.
General notes	<p>ab271915 is the carrier-free version of ab150439.</p> <p>Our carrier-free antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency.</p> <p>This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications.</p> <p>Use our conjugation kits for antibody conjugates that are ready-to-use in as little as 20 minutes with <1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.</p> <p>This product is compatible with the Maxpar[®] Antibody Labeling Kit from Fluidigm, without the need for antibody preparation. Maxpar[®] is a trademark of Fluidigm Canada Inc.</p> <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p>

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb[®] patents](#).

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Do Not Freeze.
Storage buffer	pH: 7.2 Constituent: PBS
Carrier free	Yes
Purity	Protein A purified
Clonality	Monoclonal
Clone number	EPR7412-50
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab271915 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		Use at an assay dependent concentration. Predicted molecular weight: 17 kDa.
IHC-P		Use at an assay dependent concentration.

Application notes Is unsuitable for IP.

Target

Cellular localization eIF5A2: Cytoplasm. Nucleus. Endoplasmic reticulum membrane. Nucleus > nuclear pore complex. Hypusine modification promotes the nuclear export and cytoplasmic localization and there was a dynamic shift in the localization from predominantly cytoplasmic to primarily nuclear under apoptotic inducing conditions. eIF5A: Cytoplasm. Nucleus. Endoplasmic reticulum membrane. Nucleus > nuclear pore complex. Hypusine modification promotes the nuclear export and cytoplasmic localization and there was a dynamic shift in the localization from predominantly cytoplasmic to primarily nuclear under apoptotic inducing conditions.

Images




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-eIF5A2 + eIF5A antibody [EPR7412-50] - BSA and Azide free (ab271915)

Immunohistochemical analysis of paraffin-embedded Human colon adenocarcinoma tissue labelling eIF5A2 + eIF5A with [ab150439](#) at 1/100 dilution.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

This data was developed using the same antibody clone in a different buffer formulation ([ab150439](#)).



Western blot - Anti-eIF5A2 + eIF5A antibody [EPR7412-50] - BSA and Azide free (ab271915)

All lanes : Anti-eIF5A2 + eIF5A antibody [EPR7412-50] ([ab150439](#)) at 1/1000 dilution

Lane 1 : Recombinant Human eIF5A protein ([ab87457](#))

Lane 2 : Recombinant Human eIF5A2 protein ([ab99140](#))

Lysates/proteins at 0.01 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Developed using the ECL technique.


Predicted band size: 17 kDa

Observed band size: 15,20 kDa

Exposure time: 7 seconds

Blocking buffer: 5% NFDN/TBST

This data was developed using [ab150439](#), the same antibody clone in a different buffer formulation.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-eIF5A2 + eIF5A antibody [EPR7412-50] - BSA and Azide free (ab271915)

Immunohistochemical analysis of paraffin-embedded Human testis tissue labelling eIF5A2 + eIF5A with [ab150439](#) at 1/100 dilution.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

This data was developed using the same antibody clone in a different buffer formulation ([ab150439](#)).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-eIF5A2 + eIF5A antibody [EPR7412-50] - BSA and Azide free (ab271915)

Immunohistochemical analysis of paraffin embedded Human Breast carcinoma tissue using [ab150439](#) showing +ve staining.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

This data was developed using the same antibody clone in a different buffer formulation ([ab150439](#)).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-eIF5A2 + eIF5A antibody [EPR7412-50] - BSA and Azide free (ab271915)

Immunohistochemical analysis of paraffin embedded Human Ovarian carcinoma tissue using [ab150439](#) showing +ve staining.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

This data was developed using the same antibody clone in a different buffer formulation ([ab150439](#)).



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-eIF5A2 + eIF5A antibody [EPR7412-50] - BSA and Azide free (ab271915)

Immunohistochemical analysis of paraffin embedded normal Human kidney tissue using [ab150439](#) showing +ve staining.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

This data was developed using the same antibody clone in a different buffer formulation ([ab150439](#)).

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-eIF5A2 + eIF5A antibody [EPR7412-50] - BSA and Azide free (ab271915)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors