

Product datasheet

Anti-ESGP antibody ab214934

1 Image

Overview

Product name	Anti-ESGP antibody
Description	Sheep polyclonal to ESGP
Host species	Sheep

Tested Applications & Species

Application	Species
WB	Mouse

See all applications and species data

Immunogen Recombinant full length protein corresponding to Mouse ESGP aa 1 to the C-terminus. (Derived from E.coli).
Database link: [Q2Q5T5](#)

[Run BLAST with](#)

[Run BLAST with](#)

Positive control Mouse embryonic tissue lysate.



Properties


Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Constituent: 100% PBS
Purity	Immunogen affinity purified
Purification notes	The IgG fraction of immunized serum was purified by antigen affinity chromatography.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of **ab214934** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


 Guaranteed Tested applications are guaranteed to work and covered by our Abpromise guarantee.	 Predicted Predicted to work for this combination of applications and species but not guaranteed.	 Incompatible Does not work for this combination of applications and species.
--	--	--

Application	Species
WB	 Mouse

Application	Abreviews	Notes
WB		1/1000 - 1/2000. Detects a band of approximately 12 kDa (predicted molecular weight: 10 kDa).

Target

Images



Anti-ESGP antibody (ab214934) at 1/1000 dilution + Mouse embryonic tissue lysate

Predicted band size: 10 kDa

Observed band size: 12 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors