

Product datasheet

Anti-Estrogen-induced gene 121 protein antibody ab133374

[2 Images](#)

Overview

Product name	Anti-Estrogen-induced gene 121 protein antibody
Description	Rabbit polyclonal to Estrogen-induced gene 121 protein
Host species	Rabbit
Tested applications	Suitable for: IHC-P
Species reactivity	Reacts with: Human
Immunogen	Synthetic peptide, conjugated to KLH, from internal domain of Human Estrogen-induced gene 121 protein.
Positive control	human adrenal gland, Human liver

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: 0.1% Sodium azide Constituent: 99% PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab133374** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use a concentration of 15 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

Target


Function	May play a role as a marker of hyperestrogenic state and estrogen-related type I endometrial carcinoma.
Tissue specificity	Expressed in normal endometrium but overexpressed in endometrioid tumors.
Sequence similarities	Belongs to the UPF0577 family.
Cellular localization	Membrane.

Images



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Estrogen-induced gene 121 protein antibody (ab133374)

ab133374 (15 ug/mL) was used in IHC to stain formalin-fixed, paraffin-embedded human liver, followed by biotinylated goat anti-rabbit IgG secondary antibody, alkaline phosphatase-streptavidin and chromogen.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Estrogen-induced gene 121 protein antibody (ab133374)

(ab133374 15 ug/ml) was used in IHC to stain formalin-fixed, paraffin-embedded human adrenal, followed by biotinylated goat anti-rabbit IgG secondary antibody, alkaline phosphatase-streptavidin and chromogen.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors