

## Product datasheet

# Anti-FAM63A antibody ab220932

1 Image

### Overview

<b>Product name</b>	Anti-FAM63A antibody
<b>Description</b>	Rabbit polyclonal to FAM63A
<b>Host species</b>	Rabbit

### Tested Applications & Species

Application	Species
ICC/IF	Human

See all applications and species data

**Immunogen** Recombinant fragment corresponding to Human FAM63A aa 350 to the C-terminus.  
Database link: [Q8N5J2](#)

[Run BLAST with](#)

[Run BLAST with](#)

**Positive control** A431 cells

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), 59% PBS
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of **ab220932** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**✔ Guaranteed**

Tested applications are guaranteed to work and covered by our Abpromise guarantee.

**⚠ Predicted**

Predicted to work for this combination of applications and species but not guaranteed.

**✘ Incompatible**

Does not work for this combination of applications and species.

Application	Species
ICC/IF	✔ Human
All applications	⚠ Orangutan

Application	Abreviews	Notes
ICC/IF		Use a concentration of 0.25 - 2 µg/ml. Fixation/Permeabilization: PFA/Triton X-100.

## Target


### Function

Hydrolase that can specifically remove 'Lys-48'-linked conjugated ubiquitin from proteins. Has exodeubiquitinase activity and has a preference for long polyubiquitin chains. May play a regulatory role at the level of protein turnover.

### Sequence similarities

Belongs to the peptidase MINDY family. FAM63 subfamily.

## Images



Immunofluorescent analysis of A431 cells labeling FAM63A with ab220932 at 4 µg/mL (green).

Immunocytochemistry/ Immunofluorescence - Anti-FAM63A antibody (ab220932)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors