

Product datasheet

Anti-FAM89A antibody ab188622

[3 Images](#)

Overview

Product name	Anti-FAM89A antibody
Description	Rabbit polyclonal to FAM89A
Host species	Rabbit
Tested applications	Suitable for: IHC-P
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment corresponding to Human FAM89A aa 136-162. Sequence: ALENGFFDEEEYFQEQNSLHRRDRG Database link: Q96GI7 Run BLAST with Run BLAST with
Positive control	Human colon cancer, liver and liver cancer

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.2 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab188622** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

Target


Sequence similarities Belongs to the FAM89 family.

Images




Immunohistochemical analysis of paraffin-embedded Human liver cancer labeling FAM89A with ab188622 at 1/200 dilution

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FAM89A antibody (ab188622)



Immunohistochemical analysis of paraffin-embedded Human liver labeling FAM89A with ab188622 at 1/200 dilution

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FAM89A antibody (ab188622)



Immunohistochemical analysis of paraffin-embedded Human colon labeling FAM89A with ab188622 at 1/200 dilution

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FAM89A antibody (ab188622)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors