

Product datasheet

Anti-FOXL1 antibody ab171871

1 Image

Overview

Product name	Anti-FOXL1 antibody
Description	Mouse polyclonal to FOXL1
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Recombinant full length protein corresponding to Human FOXL1 aa 1-345.

Sequence:

```

MSHLFDPRLPALAASPMLYLYGPERPGLPLAFAPAAA
LAASGRAETPQKP
PYSYALIAMA IQDAPEQRVTLNGYQFIMDRFPFYHDNR
QGWNQNSIRHN
LSLNDCFVKVPREKGRPGKGSYWTLDPRLDMFENG
NYRRRKRKPKPGPG
APEAKRPRAETHQRS AE AQPEAGSGAGGSGPAISRL
QAAPAGPSPLLDGP
SPPAPLHWPGTASPNEDAGDAAQSAAAVAVGQAAR
TGDGPGSPLRPASRS
SPKSSDKSKSFSIDSILAGKQGQKPPSGDELLGGAKP
GPGGRLGASLLAA
SSSLRPPFNASLMLDPHVQGGFYQLGIPFLSYFPLQVP
DTVLHFQ
    
```

Database link: [AAI00028.1](#)

 [Run BLAST with](#)

 [Run BLAST with](#)

Positive control	FOXL1 transfected 293T cell line lysate.
-------------------------	--

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.2 Constituent: 100% PBS

Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab171871** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 37.95 kDa.

Target

Sequence similarities Contains 1 fork-head DNA-binding domain.

Cellular localization Nucleus.

Images



All lanes : Anti-FOXL1 antibody (ab171871) at 1 µg/ml

Lane 1 : FOXL1 transfected 293T cell line lysate

Lane 2 : Non-transfected 293T cell line lysate

Secondary

All lanes : Goat Anti-Mouse IgGb at 1/2500 dilution

Predicted band size: 37.95 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors