

## Product datasheet

# GDNF Receptor alpha 1 peptide ab6099

### Overview

---

**Product name** GDNF Receptor alpha 1 peptide

### Description

---

**Nature** Synthetic

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab6099** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** Blocking

**Form** Liquid

**Additional notes** Stable for one year.

### Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Preservative: 0.02% Sodium Azide

Constituents: PBS, pH 7.2

### General Info

---

**Function** Receptor for GDNF. Mediates the GDNF-induced autophosphorylation and activation of the RET receptor.

**Sequence similarities** Belongs to the GDNFR family.

**Cellular localization** Cell membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors