

Product datasheet

Anti-GIPR antibody - C-terminal ab198694

1 References 1 Image

Overview

Product name	Anti-GIPR antibody - C-terminal
Description	Rabbit polyclonal to GIPR - C-terminal
Host species	Rabbit
Specificity	ab198694 detects endogenous level of total GIPR protein.
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Mouse, Human
Immunogen	Synthetic peptide within Human GIPR (C terminal). The exact sequence is proprietary. NCBI Accession No. NP_000155. Database link: P48546
Positive control	Mouse pancreas tissue lysate.
General notes	This product was previously labelled as Gastric Inhibitory Polypeptide Receptor

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.3 Preservative: 0.05% Sodium azide Constituents: 50% Glycerol, 49% PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab198694** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/150 - 1/1000. Predicted molecular weight: 53 kDa.

Target

Function	This is a receptor for GIP. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.
Sequence similarities	Belongs to the G-protein coupled receptor 2 family.
Cellular localization	Cell membrane.

Images



Anti-GIPR antibody - C-terminal (ab198694) at 1/150 dilution + mouse pancreas tissue lysate at 40 µg

Predicted band size: 53 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors